
# QUICK START GUIDE

## ZONA PLUS AT A GLANCE

## COMPLETING THE ZONA PLUS THERAPY

*Go to [www.zonahealth.eu](http://www.zonahealth.eu) for detailed instructions*

### *Getting Started*

Sit down in a comfortable chair and relax. Hold the device loosely in your right hand. Press the green on/off button. (Device arrives with battery pre-charged.)

### *The Calibration (Steps 1 & 2)*

1. When prompted, squeeze the device as hard as you can until "STOP" displays. Rest for 5 seconds.
2. Move the device to your left hand and follow the prompts to repeat step 1.

### *The Guided Isometric Exercise (Steps 3 - 6)*

3. With the device in your right hand, squeeze the device until the display reads "HOLD".

Continue squeezing for 120 seconds then rest for 60 seconds.

4. Move the device to your left hand, and squeeze the device until the display reads "HOLD". Continue squeezing for 120 seconds then rest for 60 seconds.
5. Move the device to your right hand and repeat step 3.
6. Move the device to your left hand and repeat step 4.

### ***The Therapy Results***

At the end of the therapy session your score is displayed. A score of 80 or above indicates a successful therapy session. To help you stay on track two pie charts are displayed showing:

How many days in the past week you used the Zona Plus successfully.

How many weeks in the past 12 you completed at least 5 therapy sessions per week.

For best results use the device a minimum of 5 days a week (unless directed otherwise by your physician). Remember it is normal to not see any results for at least 12 weeks. When your reduction in blood pressure does occur it will probably occur over a short period of time not gradually.

## **REGISTERING THE ZONA PLUS**

### **DOWNLOADING THE ZONA HEALTH APP**

1. Go to [www.zonahealth.eu/zonatherapy](http://www.zonahealth.eu/zonatherapy) and log in to your account.
2. From the "Account" tab, click "Register Zona Plus".
3. At the bottom of the dialogue box, click the "Install Zona Health Software" link and follow the setup instructions on your computer.

4. After completing step 3, using the USB cable included with the device, insert the smaller end into the bottom of the Zona Plus, and the larger into a USB port on your PC. The device will display “Charging...” and then your four character registration code.
5. Enter the characters shown on your Zona Plus screen into the field labeled “Registration Code” on the web page and click “Register”. Once registration is complete you will be prompted to unplug the device.

## UPLOADING YOUR THERAPY DATA

*Note: You must have completed the therapy at least once for the device to have data to upload.*

1. With your computer connected to the Internet, attach the device using the supplied USB cable. The Zona Plus screen should

display “Syncing...”. Once syncing is complete you may view your therapy data on the secure Zona Health web site.

## VIEWING YOUR THERAPY DATA

1. Go to [www.zonahealth.eu/zonatherapy](http://www.zonahealth.eu/zonatherapy) and log in to your account.
2. Select the “Zona Therapy” tab to view all therapy records that have been uploaded to the secure Zona Health server.

*Go to [www.zonahealth.eu/customerservice](http://www.zonahealth.eu/customerservice) for more detailed instructions.*

**WARNING:** You should not start this or any other exercise program without first consulting your doctor. To protect yourself from injury, you must read all of the warning information contained in the Zona Therapy User Guide before using the Zona Plus. Information on precautions and limitations can also be found online at: [zonahealth.de](http://zonahealth.de)

**INTENDED USE:** The Zona Plus is indicated for the treatment of hypertension.

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**Zona Health B.V.**